



Ates Medical Inc.  
50 Chemin de la Rabastalière E, Office 22A  
Saint-Bruno-de-Montarville, Québec J3V2A5, Canada  
Tél. : (450) 482-4044  
Télec. : (514) 2214079

## Hyproof Cycle Use Face Mask (Adult) GB3107



Material: 5 layers (Cotton/ Hyproof film / Polypropylene Nonwoven Fabric)  
Size: 17\*10cm


Product performance: ventilation resistance <80%, particulate matter filtration efficiency > 80%, bacteria filtration efficiency  $\geq 99.9\%$ .

Application: factories (non-special industries), kitchens, cleaning and disinfection, school institutions, outdoor activities, protection against droplets and dust.

### 5 Layers Protection

Hii-tech Hyproof Filter Layer

**BFE  $\geq 95\%$**



# Product images

8 colors



## Recommended wash methods:

1. Medical alcohol spray /  
70 degree above hot water immersion for 15 minutes.
2. Do not Machine Wash. Do not rub it.
3. Cycle Use for 200 hours.
4. Recommended wash numbers: 20 times

# Hyproof Cycle Use Face Mask (Children)



Material: 5 layers (Cotton/ Hyproof film / Polypropylene Nonwoven Fabric)  
Size: 19\*12cm (3 dimensional)


Product performance: ventilation resistance <80%, particulate matter filtration efficiency > 80%, bacteria filtration efficiency  $\geq 99.9\%$ .

Application: factories (non-special industries), kitchens, cleaning and disinfection, school institutions, outdoor activities, protection against droplets and dust.

## 5 Layers Protection

Hi-tech Hyproof Filter Layer

**BFE  $\geq 95\%$**



+ Hyproof filter +

